PNI AT198

Alcohol tester/Детектор за алкохол/Alkoholdetektor/Detector de alcohol/
Détecteur d'alcool/Alkohol detektor/Rilevatore di alcol/Alkoholdetektor/
Wykrywacz alkoholu/Detector de alcool

EN	User manual	3
BG	Ръководство за употреба	5
DE	Benutzerhandbuch	7
ES	Manual de usuario	9
FR	Manuel utilisateur	11
HU	Használati utasítás	13
IT	Manuale utente	15
NL	Handleiding	17
PL	Instrukcja obsługi	19
RO	Manual de utilizare	21

The PNI AT198 alcohol tester should only be used to indicate the possible presence of alcohol in the breath. The results obtained from the test can only be used as a basic indicator to determine whether the subject can safely drive a motor vehicle, operate equipment or perform hazardous activities.


Alcohol consumption affects a person's reflexes and judgment.

Each person responds differently to alcohol consumption, therefore the result provided by this device should only be used as a reference indicator.

Technical specifications:

- Detection range:
0.00-0.20% BAC; 0.00-2.00‰ BAC
0.00-1.00mg/L; 0.00-2.00g/L
- Accuracy: $\pm 0.02\%$ BAC- 0.05% BAC
- Preset level: 0.05% BAC; 0.50‰ BAC; 0.25mg/L; 0.50g/L
- Sensor priming time: countdown from 18 to 0
- Test duration: < 8 seconds
- LCD display with 3 digits
- Power supply: 3V DC (2 AAA alkaline batteries)
- Operating temperature: $-26^{\circ}\text{C} \sim -80^{\circ}\text{C}$
- Dimensions: 102 x 20 x 37 mm

Device introduction

1. Air intake
2. LCD screen
3. Start button
4. Loudspeaker
5. Alcohol concentration
6. The sensor primes. Wait.
7. Battery level
8. and 9. Unit for measuring the concentration of alcohol in breathed air

When the concentration of alcohol in the breathed air exceeds the preset limit, the device will emit an alarm sound.

Instructions for use

1. Press the power button for 1 second. You will hear a short confirmation sound.
2. The screen will display the message WAIT and the countdown from 18 to 0, during which the sensor primes.

Note: If the device has not been used for a long time, the sensor priming time may be longer.

3. After the sensor has primed, the WAIT message disappears from the screen and the “C” and “blow” icons will appear. From this moment you can start testing.
4. Take a deep breath and blow for 3-4 seconds near the holes above the screen. Keep blowing until you hear a sound. The 0.00 icon will flash on the screen, and the “blow” message will disappear.
5. In a short time, the test result will be displayed on the screen accompanied by a sound signal. The result remains on the screen for a few seconds, after which the device will close automatically.
6. If you want to do another test, press the start button again and repeat the steps above.

Warnings

- Do not clean the device with chemical cleaning agents, as they release gases that could affect the performance of the sensor.
- If the device has not been used for a long time, it is possible that the first test will not succeed. In this case, repeat the test procedure.
- To ensure a correct result, take the test at least 10 minutes after consuming alcohol.
- Do not use the device in environments with aggressive gases, such as chlorine.
- Under normal conditions of use, the lifetime of the sensor is approximately 2 years.

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