## PNI FG8118

Muscle massage gun / Устройство за вибромасаж /  
Vibrationsmassagerät / Dispositivo de vibromasaje / Appareil  
de vibromassage / Vibromassázs készülék / Dispositivo di  
vibromassaggio / Vibromassage apparaat / Urządzenie do masażu  
wibracyjnego / Aparat de vibromasaj

EN	User manual .....	3
BG	Ръководство за употреба .....	5
DE	Benutzerhandbuch .....	7
ES	Manual de usuario .....	9
FR	Manuel utilisateur .....	11
HU	Használati utasítás .....	13
IT	Manuale utente .....	15
NL	Handleiding .....	17
PL	Instrukcja obsługi .....	19
RO	Manual de utilizare .....	21

## Battery charging

Before the first use, please charge the battery with the USB power cable included in the package. Do not use the massage gun while it is charging.

## Turn on/off

Press the ON/OFF button to turn on the massage gun and press it again to turn it off.

## Safety warnings

- Clean the massage gun with a slightly damp cloth, then wipe it well with a dry cloth.
- Keep the massage gun in well-ventilated space, away from excessive heat and humidity.
- In case of malfunction, do not try to repair the massage gun yourself. Call a specialized service center.
- After charging the battery, remove the power adapter.
- Do not leave the massage gun within the reach of children.
- Do not submerge the massage gun in water.
- If you feel pain or discomfort, stop using this device immediately.
- Use this massage gun only on soft areas of the body. Do not try to massage the head area or any other area of the body

with visible bones.

- Use the massage gun by gently pressing and moving it on the surface of the body, but not more than 60 seconds on one area.
- Do not use this massage gun if you have wounds or bruises on your body.
- To charge the massage gun, use only the power adapter included in the package.
- Use only the massage accessories included in the package.
- If you suffer from serious diseases, consult your doctor before using this device.

## Technical specifications

- Power supply: Li-Ion 1200 mAh 7.4V battery
- Supply voltage: 5V
- Speed: 6 levels
- Power: 20W
- Battery level indicator
- Dimensions: 65 x 145 x 160 mm
- Working temperature: -26°C ~ +80°C

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