PNI FG8121

Muscle massage gun / Устройство за вибромасаж /
Vibrationsmassagerät / Dispositivo de vibromasaje / Appareil
de vibromassage / Vibromassázs készülék / Dispositivo di
vibromassaggio / Vibromassage apparaat / Urządzenie do masażu
wibracyjnego / Aparat de vibromasaj

EN	User manual	3
BG	Ръководство за употреба	5
DE	Benutzerhandbuch	7
ES	Manual de usuario	9
FR	Manuel utilisateur	11
HU	Használati utasítás	13
IT	Manuale utente	15
NL	Handleiding	17
PL	Instrukcja obsługi	19
RO	Manual de utilizare	21

Battery charging

Before the first use, please charge the battery with the power adapter included in the package. Do not use the massage gun while it is charging. The charge level of the battery will be displayed on the screen at the bottom of the massage gun.

Turn on/off

Move the ON/OFF button to the ON position to turn on the massage gun and to the OFF position to turn it off.

Safety warnings

- Clean the massage gun with a slightly damp cloth, then wipe it well with a dry cloth.
- Keep the massage gun in well-ventilated space, away from excessive heat and humidity.
- In case of malfunction, do not try to repair the massage gun yourself. Call a specialized service center.
- After charging the battery, remove the power adapter.
- Do not leave the massage gun within the reach of children.
- Do not submerge the massage gun in water.
- If you feel pain or discomfort, stop using this device immediately.
- Use this massage gun only on soft areas of the body. Do not

try to massage the head area or any other area of the body with visible bones.

- Use the massage gun by gently pressing and moving it on the surface of the body, but not more than 60 seconds on one area.
- Do not use this massage gun if you have wounds or bruises on your body.
- To charge the massage gun, use only the power adapter included in the package.
- Use only the massage accessories included in the package.
- If you suffer from serious diseases, consult your doctor before using this device.

Technical specifications

- Power supply: Li-Ion 1500 mAh battery
- Supply voltage: 12V DC
- Speed: 12 levels
- Power: 20W
- Battery level indicator
- Dimensions: 175 x 220 x 65 mm
- Working temperature: -26°C ~ +80°C

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